## VINTAGE HALO MARGIN PORCELAIN

Margin porcelain for metal bond

VINTAGE HALO MARGIN PORCELAIN can reproduce an esthetic appearance close to the natural teeth easily and quickly.

### [NOTE ON USE]

1. While working on the finishing and polishing, vacuum cleaner or dustproof mask is recommended not to inhale the dust.
2. For finishing and polishing, eye protection glass is recommended.
3. This product is to use by dental professional only.
4. The product shall not be used except for the intended use.

### [INSTRUCTIONS ON USE]

1. Preparation of abutment tooth with shoulder margin.
2. Metal free margin at shoulder part.
3. After opaque firing, take out the crown from the die and apply Halo Margin porcelain with slight condensation, then replace the crown to the die firmly.
4. Compress the moist powders to the shoulder with carver.
5. Make sure that crown can be removed from the die properly and fit to the die again.
6. Smooth the surface with dry brush to eliminate the excess powder, then proceed to firing.
7. After first firing, remove the crown from the die, and apply small amount of Margin porcelain additionally to the shoulder. It would be easier to build up when using the dilution of CPM liquid with distilled water at 1:1 ratio.
8. Lightly compress the crown to the die to extrude the excess powder, which is to be taken away by brush. Then twist the crown gently to take off from the die, and proceed to the second firing.
9. After second firing, build up Vintage Halo body and incisal powders, and contour the shape.
10. After biscuit bake, proceed to contour adjustment and self-glazing.

### ■ Should the Margin require an addition after glazing, use CPM or CPM fine powder for correction.

CPM...If fitting condition was fairly poor.

CPM fine...If fitting condition is not good.

1. Remove the crown from the die, and apply small amount of CPM or CPM fine powder mixed with VINTAGE CPM Modelling Liquid to the tip of margin, then slightly condense.
2. Lightly compress the crown to the die to extrude the excess powder, which is to be taken away by brush.
3. Twist the crown gently to take off the die, and proceed to the firing according to the schedule.

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Correction powder (2 colors) : CPM, CPM FINE

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
**【築盛図・Layering structure・Schichtsystem】**

**【TEMPERATURE RISE RATE: 50 – 60°C/min.】**

(NOTE) Trial firing is recommended before substantial use to obtain the precise temperature, as the temperature will vary depending on the type of furnace. The recommended firing schedule refers to porcelain furnaces with an operational voltage of 110V. Porcelain furnaces with an operational voltage of 220-240V can cause 10-20°C lower final temperature at firing 3, 4 and 5.

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Manufactured by  
**SHOFU INC.** 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan  
 SHOFU INC. Singapore Branch 10 Science Park Road #03-12, The Alpha, Science Park II, Singapore 117684  
 SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA  
 SHOFU DENTAL PRODUCTS LTD. Duke's Factory, Chiddingstone Causeway, Tonbridge, Kent, TN11 8JU, UK  
 SHOFU DENTAL SUPPLIES (Shanghai) CO. LTD. No.645 Jiye Road, Sheshan Industry Park, Songjiang 201602, Shanghai, China  
 [EC][REP] SHOFU DENTAL GmbH Am Brühl 17, 40878 Ratingen, Germany www.shofu.de