



Ceravety Press & Cast

Universal Speed Investment for Pressable Ceramics and Cast Alloys

Directions for Use

FEATURES

Ceravety Press & Cast is a universal speed investment for pressing / press-over and casting procedures, but also suitable for conventional overnight heating. It provides a long working time of 5-6 minutes. Thanks to its excellent expansion control, Ceravety Press & Cast optimally compensates the shrinkage of non-precious and precious alloys and makes the surfaces of pressed and casted objects very smooth. Ceravety Press & Cast features easy divestment, minimizing the stress for pressed objects with low wall thicknesses. When using the press-over technique, metal or ceramic objects can also be invested.

This product has been specifically optimized for the use of IPS e.max Press lithium disilicate ceramic (Ivoclar AG); it produces only very thin reaction layers.

The recommended **Powder : Liquid mixing ratio is 100 g : 20 ml**. The expansion of the investment can be precisely controlled by changing the Liquid concentration with distilled water.

INDICATIONS

- Suitable for all commercially available pressing pellets (high and low melting)
 - IPS e.max Press lithium disilicate ceramic (Ivoclar AG)
 - Non-precious and precious alloys
-

NOTES FOR DENTAL PERSONNEL

- This product contains free silica; long-term inhalation of free silica may cause lung damage. Use local dust extractor, dust protection mask etc. to avoid harmful influence of the dust on the human body.
 - Use this product only in well-ventilated rooms, equipped with a suitable ventilation system or a fan, to avoid inhalation of the gas produced by heating the material.
 - When grinding this product, use protective equipment, such as protective glasses.
 - Avoid contact with eyes. In case of contact with eyes, immediately flush eyes with plenty of water and seek medical assistance.
 - Use protective gloves to avoid direct contact with this product, because both the Liquid and the Powder are alkaline materials.
 - This product is intended for use by dental professionals only.
-

SPEED TECHNIQUE / PRESSABLE CERAMICS

- Wax-mount the wax pattern on the base of the pressing cylinder (muffle) in the usual way and add a paper or silicone ring.
- When using the press-over technique, metal or ceramic objects can also be invested.
- Mix the Powder and the Liquid in the correct ratio, as described below.
- Carefully fill the muffle with this mixture, using the lowest vibration setting.

Working time: 5-6 min (at 23°C room temperature)

Mixing time: Hand mixing for 30 sec, then vacuum mixing for 60 sec. Pressure investment is not necessary.

Setting time: 20 to max. 25 min, including the mixing procedure. Remove the base and the muffle former 15 min after mixing and allow the muffle to dry. Then place it in a furnace preheated to the final temperature.

Preheating temperature: 850°C

Pressable ceramics	Indications	Liquid : distilled water 100 g muffle	Liquid : distilled water 200 g muffle
Mixing ratio	Inlays / Onlays	5 ml : 15 ml	10 ml : 30 ml
	Crowns / press-over technique	8 ml : 12 ml	16 ml : 24 ml
	Bridges	8 ml : 12 ml	16 ml : 24 ml
Holding time at final temperature		45 min	60 min

Divesting: After cooling, remove the investment from the pressed object, using a cutting disk or glass beads (50 µm at approx. 4 bar), without exposing the objects. Then reduce the pressure to 1-2 bar and carefully expose the objects.

Note: When pressing lithium disilicate, a thin reaction layer may form, depending on the furnace used. This layer must be removed by sandblasting with aluminum oxide (50-100 µm at approx. 1-2 bar); please observe the manufacturer's instructions!

CONVENTIONAL PREHEATING / PRESSABLE CERAMICS

Preheating times: After setting, place the muffle in a cold furnace. The heating rate should be approx. 3-5°C per minute in a furnace with linear control.

Note: When preheating overnight, the investment must be wrapped in cellophane to prevent it from drying out.

Preheating times	100 g muffle	200 g muffle	300 g - 480 g muffle
Holding time at 260°C	30 min	45 min	70 min
Holding time at 580°C	20 min	30 min	60 min
Holding time at final temperature	40 min	45 min	70 min

SPEED TECHNIQUE / CAST ALLOYS

- Wax-mount the wax pattern on the muffle base former in the usual way and add the muffle ring.
- Mix the Powder and the Liquid in the correct ratio, as described below.
- Carefully fill the muffle with this mixture, using the lowest vibration setting.

Working time: 5-6 min (at 23°C room temperature)

Mixing time: Hand mixing for 30 sec, then vacuum mixing for 60 sec. Pressure investment is not necessary.

Setting time: 20 to max. 25 min, including the mixing procedure. Then place the muffle in a furnace preheated to the final temperature.

Preheating temperature: 850°C to max. 900°C (according to the alloy manufacturer's instructions)

Note: These values apply to a **Powder : Liquid mixing ratio of 100 g : 20 ml**.

Alloys	Indications	Liquid : distilled water 1 x muffle	Liquid : distilled water 3 x muffle
Non-precious alloy	Crowns	16 ml : 4 ml	32 ml : 8 ml
	Bridges	18 ml : 2 ml	36 ml : 4 ml
Precious alloy	Crowns	10 ml : 10 ml	20 ml : 20 ml
	Bridges	12 ml : 8 ml	24 ml : 16 ml
Holding time at final temperature		45 min	60 min

CONVENTIONAL PREHEATING / CAST ALLOYS

Preheating times: After setting, place the muffle in a cold furnace. The heating rate should be approx. 3-5°C per minute in a furnace with linear control.

Note: When preheating overnight, the investment must be wrapped in cellophane to prevent it from drying out.

Preheating times	1 x muffle	3 x muffle	6 x - 9 x muffle
Holding time at 260°C	30 min	45 min	70 min
Holding time at 580°C	20 min	30 min	60 min
Holding time at final temperature	40 min	45 min	70 min

PHYSICAL PROPERTIES

Basic Data

Powder : Liquid mixing ratio	100 g : 20 ml
Working time	6 min
Setting time	9 min
Thermal expansion / 850°C	0.9 %
Compressive strength after 2 hours	10.0 MPa

Adjustment of thermal expansion by means of Liquid concentration

Liquid concentration (%)	Setting expansion	Thermal expansion	Total expansion
0			
20	0.40 %	0.60 %	1.00 %
40	0.50 %	0.75 %	1.25 %
60	0.60 %	0.80 %	1.40 %
80	0.70 %	0.85 %	1.55 %
100	1.00 %	0.90 %	1.90 %

COMPOSITION

POWDER: Silica, quartz, cristobalite powder and others

LIQUID: Colloidal silica and water

STORAGE

Ceravety Press & Cast should be stored in a dry place at room temperature. Do not store the Powder at high humidity and the Liquid in direct sunlight. The Liquid is frost-sensitive and must not be stored at temperatures below 0°C!

PACKAGING

Ceravety Press & Cast Powder (PN 6966): 120 x 100 g

Ceravety Press & Cast Powder (PN 6968): 30 x 100 g

Ceravety Press & Cast Liquid (PN 6967): 2 liters

Ceravety Press & Cast Liquid (PN 6969): 300 ml



Manufacturer

SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan

SHOFU DENTAL ASIA-PACIFIC PTE. LTD. 10 Science Park Road, #03-12, The Alpha, Science Park II, Singapore 117684

SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, California 92078-4059, USA

SHOFU UK Riverside House, River Lawn Road, Tonbridge, Kent, TN9 1EP, UK

SHOFU DENTAL TRADING (SHANGHAI) CO., LTD. No. 645 Jiye Road, Sheshan Industrial Park, Songjiang, Shanghai 201602, China

EC REP SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany (www.shofu.de)