RESIN GLAZE

EN Instructions for use	DE Gebrauchsanleitung	FR Mode d'emploi	NL Gebruiksaanwijzing
ES Instrucciones de uso	IT Istruzioni per l'uso	RU Инструкция по применению	PL Instrukcja obsługi
RO Instrucțiuni de utilizare	SV Bruksanvisning	EL Οδηγίες χρήσης	PT Instruções de uso

02:2019-03
71599

	Manufacturer SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan www.shofu.com	EC REP SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, California 92078-4059, USA	SHOFU DENTAL ASIA-PACIFIC PTE.LTD. 10 Science Park Road, #03-12, The Alpha, Science Park II, Singapore 117684	

EN READ CAREFULLY BEFORE USE

LIGHT-CURED SURFACE GLOSS COATING MATERIAL

PRODUCT DESCRIPTION

SHOFU RESIN GLAZE is a light-cured, temporary surface gloss coating material for restorations. Liquid, glazing material can be bonded more securely when used in conjunction with Primer.

INDICATIONS FOR USE

SHOFU RESIN GLAZE is indicated to temporarily glaze the following restoration surfaces:

- Artificial resin teeth
- Acrylic denture base materials
- CAD/CAM PMMA materials
- CAD/CAM hybrid ceramic materials
- Indirect composite resins
- Provisional materials

PRECAUTIONS AND WARNINGS

- Do not use this product on patients with known allergies to this material and/or methacrylate monomer.
- Operators with known allergies to this material and/or methacrylate monomer must not use this product.
- If any inflammation or other allergic reactions occur on either the patient or operator, immediately discontinue use and seek medical advice.
- Avoid contact with skin and eyes. In case of accidental contact, immediately flush with plenty of water and seek medical advice.
- Since this material is not permanently durable, use it as temporary glazing material.
- When applying this material to the clasp of flexible removable partial denture, cracks or removal might occur when placing or removing the denture.
- Do not use this product for any purposes other than specifically outlined in these instructions for use.
- Use this product within the expiration date indicated on the package. (Example ☞ YYYY-MM-DD → Year-Month-Date of the expiration date)
- This product is intended for use by dental professionals only.

DIRECTIONS FOR USE

USE EXTRAORALLY

- Preparation of the surface to be glazed
Contour and finish the surface of the restoration to be glazed with an appropriate abrasive (diamond point, silicone point, carborundum point or carbide bur etc.). Sandblast with 50 µm aluminum oxide (0.1-0.2 MPa/1-2 bar) depending on the case. Then, clean with alcohol and dry.
- Application
SHOFU RESIN GLAZE Primer can be applied before Liquid (glazing material) to enhance bonding to the restoration.
 - Application of Primer
Apply Primer on the prepared surface with a brush, then dry for 10 seconds.
 - Application of Liquid (glazing material)
Apply Liquid to the treatment surface in a thin layer with a brush. (Use separate brushes for Primer and Liquid.) Don't air blow.
- Light-curing
Light-cure for 3 minutes with a dental lab curing unit with Halogen lamp.
 - Curing units having only Blue LED light source (wavelength of between 450 nm and 470 nm) cannot be used for this material.
 - Applicable wavelength for this material is between 400 nm and 410 nm regardless of light source.

NOTE ON USE

- Use this material in dust-free area.
- Do not mix this material with other materials.
- When storing these materials in a refrigerator, allow them to reach room temperature before use.
- Thoroughly remove resin dust, oil or other contaminants from the resin surface before application.
- Always use neat brush and move it gently in one direction to apply this material. Frayed or worn brush causes entrapment of air bubbles.
- Avoid applying too thickly to avoid entrapping air bubbles or degrading the physical properties.
- Avoid exposure to intense light, such as direct sunlight or work light etc. to avoid unintended polymerization of the material.
- Clean the brush with ethanol after applying this material.
- If leaving the liquid uncapped, volatile agent might evaporate, causing solidification or contamination of the material. Cap the bottle securely immediately after each use.
- Refer to dental lab light-curing unit manufacturer's instructions regarding their specific wavelength and light intensity to light-cure this material properly.
- Unless this material is light-cured sufficiently, uncured layer causes staining or discoloration. Light-cure for the time duration recommended by the manufacturer.
- Follow the instructions for use of each dental material, instrument or equipment to be used in conjunction with this product.
- In order to ensure proper polymerization of this material, exercise caution to find any changes in light intensity, dirt or lamp life.

COMPOSITION

Liquid: Multi-functional monomer, Methacrylate monomer and others
Primer: Multi-functional monomer, Methacrylate monomer and others

STORAGE

- Store in dark place at room temperature (1-25 °C/34-77 °F). Avoid high temperature and humidity. Keep away from direct sunlight and any source of ignition.
- Tightly close the cap after each use.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

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