Universal Opaque

EN Instructions for use
ES Instrucciones de uso
RO Instrucțiuni de utilizare

DE Gebrauchsanleitung
IT Istruzioni per l'uso
SV Bruksanvisning

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SHOFU DENTAL GmbH
Am Brühl 17, 40878 Ratingen, Germany

SHOFU DENTAL CORPORATION
1225 Stone Drive, San Marcos,
California 92078-4059, USA

SHOFU DENTAL ASIA-PACIFIC PTE. LTD.
10 Science Park Road, #03-12,
The Alpha, Science Park II, Singapore 117684

03:2019-03
71618

EN Light-curing opaque paste

READ CAREFULLY
BEFORE USE

PRODUCT DESCRIPTION

SHOFU Universal Opaque consists of two articles, PRE-OPAQUE and OPAQUE, to be used in conjunction with SOLIDEX, CERAMAGE and other crown and bridge composite resins.

INDICATIONS FOR USE

- Masking the metal color of facing crown
- Masking the substructure color of jacket crown and laminate veneer etc.
- Masking the substructure color of restorations and prostheses made of porcelain or crown and bridge composite resin etc. when repairing them intraorally or extraorally

PRECAUTIONS AND WARNINGS

- Do not use this product on patients with known allergies to this product and/or methacrylate monomer.
- Operators with known allergies to this product and/or methacrylate monomer must not use this product.
- If any inflammation or other allergic reactions occur on either the patient or operator, immediately discontinue use and seek medical advice.
- Use protective plastic medical (dental) gloves and glasses to avoid sensitization to this product. In case of accidental contact with oral soft tissue or skin, immediately blot with alcohol moistened cotton ball, and rinse with plenty of water. In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
- There have been some reports of operators developing sensitization reaction by handling uncured material that incorporates some chemical components included in this product with bare hands.
- Follow the instructions for use of each dental material, instrument or equipment to be used in conjunction with this product.
- Do not use this product for any purposes other than specifically outlined in these instructions for use.
- Use this product within the expiration date indicated on the package. (Example ☞ YYYY-MM-DD → Year-Month-Date of the expiration date)
- This product is intended for use by dental professionals only.

DIRECTIONS FOR USE

SHOFU Universal Opaque can be used for a variety of cases to mask the substructure. Thus, fabrication of resin-based facing crown and bridge is shown as a typical indication of this product.

- Preparation of frameworks
 - Produce a metal frame in the conventional manner. Sandblast the metal surface with alumina sand. Then, apply metal primer (M.L. Primer, SHOFU Universal Primer etc.) to the surface.
- Application of PRE-OPAQUE
 - 2-1. Dispense adequate amount of PRE-OPAQUE onto a paper pad, dish or directly onto a brush.
 - 2-2. Apply the dispensed material to the pretreated metal surface.
 - 2-3. Light-cure the applied PRE-OPAQUE using a recommended laboratory light-curing unit for the prescribed time.
- Application of OPAQUE
 - 3-1. Dispense adequate amount of OPAQUE onto a paper pad, dish or directly onto a brush.
 - 3-2. Apply the dispensed material to the PRE-OPAQUE applied surface.
 - 3-3. Light-cure the applied OPAQUE using a recommended laboratory light-curing unit for the prescribed time.

- When using a laboratory light-curing unit other than Solidilite V, ensure optimal light intensity to light-cure these materials properly in reference to the instructions of each light-curing unit.
- Repeat the application and light-cure process if necessary.
- If applying to a large area or multiple teeth, divide the area into some parts to light-cure and change the direction of light irradiation.
- If applying in multiple layers, light-cure each layer properly.
- After confirming that OPAQUE is light-cured completely, build up crown and bridge composite resin onto the cured OPAQUE, complete the crown form and light-cure the composite. Contour, finish and polish the restoration in the conventional manner.

NOTE ON USE

- Tightly close the cap immediately after each use.
- Use an appropriate brush for application. Clean the used brush with alcohol etc.
- Shield the dispensed material from ambient light using a light shielding cover.
- Avoid handling SHOFU Universal Opaque in the area exposed to intense light such as close to the window or immediately beneath the work light to avoid unintended polymerization of the paste.
- Use separate brushes for PRE-OPAQUE and OPAQUE.

COMPOSITION

PRE-OPAQUE: UDMA, Silica powder, Pigments, Polymerization initiator and others
OPAQUE: UDMA, Silica powder, Pigments, Polymerization initiator and others

STORAGE

Store in dark place at room temperature (1-25 °C). Avoid high temperature and humidity. Keep away from direct sunlight and any source of ignition.

SHADES

PRE-OPAQUE (1 shade)	PRE-OPAQUE
OPAQUE (31 shades)	A10, A20, A30, A3.50, A40, B10, B20, B30, B40, C10, C20, C30, C40, D20, D30, D40, rootAO, R20, R30, R3.50, W0, W20, W30, GUM-O, BGO, GO, MO, InO, VpO, VrO, PO

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

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