

# Perfect Gloss for Porcelain and Composite



## CeraMaster® • CompoMaster®

The unique blend of silicone and the even distribution of diamond particles as well as a variety of popular shapes for polishing porcelain and composite surfaces give restorations a superb shine.

**CeraMaster Coarse** and **CompoMaster Coarse** feature the same positive properties as the proven **CeraMaster** and **CompoMaster** polishers, but due to their bigger particle size these instruments are suitable for a more efficient contouring and pre-polishing.

Thanks to their shape durability and long tool life these polishers are economical; moreover they eliminate the need for expensive polishing pastes.



Proven Products for Better Dentistry

# The diamond impregnated professional polishers

## CeraMaster® Coarse • CeraMaster®

The **CeraMaster Coarse** polishers are designed for a quick finishing and pre-polishing of porcelain restorations and provide a smooth surface.

They particularly show their forte when adjusting ceramic coping margins made from materials such as aluminium oxide and zircon oxide.

For superpolishing **CeraMaster** polishers are recommended.



ISO No.	Shape	Bullet	Knife	Cup	Mini-Point	KN7	WH6	unmounted Wheel*
		040	090	065	030	125	120	220
CeraMaster Coarse	CA 3 pcs	121C	122C	123C	124C			
	HP 3 pcs	126C			129C	127C	128C	132C
CeraMaster	CA 3 pcs	0121	0122	0123	0124			
	HP 3 pcs	0126			0129	0127	0128	

\* 1 piece + mandrel / box

### Additionally assorted kits are available:

**CeraMaster Coarse** – assorted kit, CA shank: 1x Bullet, Knife, Cup, Mini-Point

PN 125C

**CeraMaster Coarse** – assorted kit, HP shank: 1x Bullet, Mini-Point, KN7, WH6

PN 130C

**CeraMaster** – assorted kit, CA shank: 1x Bullet, Knife, Cup, Mini-Point

PN 0125

**CeraMaster** – assorted kit, HP shank: 1x Bullet, Mini-Point, KN7, WH6

PN 0130

## CompoMaster® Coarse • CompoMaster®

**CompoMaster Coarse** polishers have been designed for contouring and prepolishing rough composite surfaces in one step.

Finally the prepared surfaces can be superpolished with the **CompoMaster** instruments.



ISO No.	Shape	Bullet	Knife	Cup	Mini-Point	KN7	WH6	unmounted Wheel*
		040	090	065	030	125	120	220
CompoMaster Coarse	CA 3 pcs	111C	112C	113C	114C			
	HP 3 pcs	116C			119C	117C	118C	131C
CompoMaster	CA 3 pcs	0111	0112	0113	0114			
	HP 3 pcs	0116			0119	0117	0118	

\* 1 piece + mandrel / box

### Additionally assorted kits are available:

**CompoMaster Coarse** – assorted kit, CA shank: 1x Bullet, Knife, Cup, Mini-Point

PN 115C

**CompoMaster Coarse** – assorted kit, HP shank: 1x Bullet, Mini-Point, KN7, WH6

PN 120C

**CompoMaster** – assorted kit, CA shank: 1x Bullet, Knife, Cup, Mini-Point

PN 0115

**CompoMaster** – assorted kit, HP shank: 1x Bullet, Mini-Point, KN7, WH6

PN 0120

**SHOFU DENTAL GMBH**

Am Brüll 17 · D-40878 Ratingen · Germany  
Phone +49 (0) 2102 / 86 64-0 · Fax +49 (0) 2102 / 86 64-64  
E-mail: info@shofu.de · www.shofu.de

