


# Schicht-System

- Opal Inzisal
- Body
- Opaque-Dentin
- Opaque
- Metall



## Farbtabelle

Farbe		A1	A2	A3	A3.5	A4	root A
Opaque		A1O	A2O	A3O	A3.5O	A4O	root AO
Opaque-Dentin		OD-A1	OD-A2	OD-A3	OD-A3.5	OD-A4	OD-root A
Body		A1 B	A2 B	A3 B	A3.5 B	A4 B	root A B
Opal Inzisal oder Standard Inzisal	57	○					
	58	■	○ ■				
	59			○ ■	○ ■		
	60					○ ■	○ ■
Transluzent	T	■	■	■	■	■	■

○ Opal Inzisal


■ Standard Inzisal

B1	B2	B3	B4	root B
B1 O	B2 O	B3 O	B4 O	root B O
OD-B1	OD-B2	OD-B3	OD-B4	OD-root B
B1 B	B2 B	B3 B	B4 b	root B B
○				
	○ ■			
		○ ■		
			○ ■	○ ■
■	■	■	■	■

C1	C2	C3	C4	rootC	D2	D3	D4
C1 O	C2 O	C3 O	C4 O	root C O	D2 O	D3 O	D4 O
OD-C1	OD-C2	OD-C3	OD-C4	OD-root C	OD-D2	OD-D3	OD-D4
C1 B	C2 B	C3 B	C4 B	root C B	D2 B	D3 B	D4 B
○ ■	○ ■				○		
		○ ■			■	○	○ ■
			○ ■	○ ■		■	
■	■	■	■	■	■	■	■

# Build-up System

- Opal Incisal
- Body
- Opaque-Dentine
- Opaque
- Metal



# Color-Chart

Shade		A1	A2	A3	A3.5	A4	root A
Opaque		A1O	A2O	A3O	A3.5O	A4O	root AO
Opaque-Dentine		OD-A1	OD-A2	OD-A3	OD-A3.5	OD-A4	OD-root A
Body		A1 B	A2 B	A3 B	A3.5 B	A4 B	root A B
Opal Incisal or Standard Incisal	57	○					
	58	■	○ ■				
	59			○ ■	○ ■		
	60					○ ■	○ ■
Translucent	T	■	■	■	■	■	■

○ Opal Incisal

■ Standard Incisal

B1	B2	B3	B4	root B
B1 O	B2 O	B3 O	B4 O	root B O
OD-B1	OD-B2	OD-B3	OD-B4	OD-root B
B1 B	B2 B	B3 B	B4 b	root B B
○				
	○ ■			
		○ ■		
			○ ■	○ ■
■	■	■	■	■

C1	C2	C3	C4	rootC	D2	D3	D4
C1 O	C2 O	C3 O	C4 O	root C O	D2 O	D3 O	D4 O
OD-C1	OD-C2	OD-C3	OD-C4	OD-root C	OD-D2	OD-D3	OD-D4
C1 B	C2 B	C3 B	C4 B	root C B	D2 B	D3 B	D4 B
○ ■	○ ■				○		
		○ ■			■	○	○ ■
			○ ■	○ ■		■	
■	■	■	■	■	■	■	■