

MZ

Primer Plus

Strong Bonds without Mechanical Retention



SHOFU MZ Primer Plus

Easy, Efficient and Universal

The new MZ Primer Plus creates extremely reliable bonds between resins and all dental alloys, zirconia or alumina.

Without the use of retention beads, the unique primer ensures durable bond strengths when used in combination with the light-cured opaque pastes of the SHOFU Universal Opaque System, even in a moist environment.

Indications

- Application of light-cured opaque resins to frameworks
- Manufacturing and repair of C&B restorations
- Intraoral and extraoral repair of metal-ceramic restorations
- Bonding of metal bases to denture acrylic resin
- Bonding of resin to orthodontic appliances



SHOFU MZ Primer Plus




SHOFU Universal Pre-Opaque

Benefits

- Easy and quick handling without any additional equipment
- Transparent, no shadows or dark margins in the bonding area
- Reliable bonds – even without the use of retention beads
- For all dental alloys and zirconia-based frameworks
- For intraoral and extraoral use

SHOFU MZ Primer Plus, the successor to the clinically proven M.L. Primer, may be used for a wide range of substrates, including all precious and nonprecious alloys, zirconia and alumina, and is characterised by easy and efficient handling without any further accessories or heat treatments. It is the universally reliable solution for bonding resins to all frameworks made by CAD/CAM technology.



SHOFU Universal Opaque

The combination of SHOFU Universal Opaque and MZ Primer Plus offers a coordinated bonding system. These bondings are extremely resistant to hydrolysis and thermal loadings.

Order Information



SHOFU MZ Primer Plus

Contents 5 ml, PN Y0006

Related Products



SHOFU Universal Opaque

Light-curing paste opaques



Ceramage

Light-curing micro-hybrid composite



Ceramage UP

Flowable, light-curing micro-hybrid composite for freehand layering technique



SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan (www.shofu.co.jp)
SHOFU DENTAL ASIA-PACIFIC PTE. LTD. 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)
SHOFU DENTAL TRADING (SHANGHAI) CO., LTD. No. 645 Jiye Road, Sheshan Industrial Park, Songjiang, Shanghai 201602, China (www.shofu.com.cn)
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)
SHOFU UK Riverside House, River Lawn Road, Tonbridge, Kent, TN9 1EP, UK (www.shofu.co.uk)
SHOFU DENTAL GmbH [EC REP] Am Brühl 17, 40878 Ratingen, Germany (www.shofu.de)

SHOFU is a registered trademark of SHOFU INC. All other trademarks and registered trademarks are the property of their respective holders. SHOFU INC. reserves the right to change specifications without notice.