

All Ceramic Preparation Kit



Preparation system for all-ceramic restorations

The assortment of the 17 instruments is ideally suited for the preparation of the anterior area, for the preparation of premolars and molars as well as the inlay and facet preparation. The shapes of the abrasives are optimally designed for the aspired preparation, ensuring a shape congruent performance.

Contents

Robot Points - 1 each:

808V-1, 0825-1, 825V-1, 0826-1, 826V-1, 0833-1,
833V-1, 0836-1, 836V-1, 0839-1, 839V-1, 0893-1,
893V-1, 0896-1, 0897-1, 0899-1, 899V-1

PN	Quantity	Product
0955	1 pc.	All Ceramic Preparation Kit

