

Amalgam Polishing Kit

Polishing kits for amalgam and precious alloys

The kits are available in contra angle or friction grip shank and contain all important instruments for all polishing steps up to high gloss.

Contents

Cup: 3x Brownie, 2x Greenie, 1x Supergreenie

Minipoint: 3x Brownie, 2x Greenie, 1x Supergreenie

PN	Quantity	Product
0308	1 pc.	Amalgam Polishing Kit, CA shank
0309	1 pc.	Amalgam Polishing Kit, FG shank

