




BEAUTIFIL Flow Plus

Injectable anterior and posterior filling composite

Beautifil Flow Plus combines the application properties of a flow with the stability, durability and aesthetic of packable composites. It is approved for every classification – including occlusal surfaces and proximal margins. With improved handling properties allowing an injectable application, creating posterior restorations has never been easier. Moreover, the ability of Beautifil Flow Plus to release and recharge fluoride and its resistance to bacterial colonisation makes it ideal for the prevention of secondary caries.



Beautifil Flow Plus is available in two viscosities, "F00" is extremely stable and "F03" is moderately flowable.

Both are suited to the known flow indications, but as well to the occlusion-bearing areas which means they can also be used for class I and II fillings. With the 360° rotary finger grip and a replaceable fine cannula the Flow Plus composites are easily applied directly out of the syringes.

A wide range of colours and various opacities are available for aesthetic restorations in single and multi shade technique.

BENEFITS

- Approved for restorations of all cavity classes
- High tensile and flexural strength
- "Cone technique" application
- Natural aesthetics with chameleon effect
- Sustainable fluoride release and recharge
- Exceptional radiopacity

Indications

- F00 – ZERO FLOW
Shaping of occlusal anatomy, marginal ridges and intricate surface details, restorations of all cavity classes
- F03 – LOW FLOW
Restoration of pits, gingival boxes for surface modification; serves as a baseline





F00 – absolutely stable (Zero Flow)

Shade	PN 2.2 g	Shade	PN 0.21 g
SYRINGE		TIPS, 20 pcs.	
A1	2001	A1	2201
A2	2002	A2	2202
A3	2003	A3	2203
A3.5	2004	A3.5	2204
A4	2005	A4	2205
B1	2062	B1	2206
B2	2063	B2	2207
C2	2064	C2	2208
D2	2065	D2	2209
SA0.5	2000	SA0.5	2200
SA0.50 (opacious A0.5)	2006	SA0.50 (opacious A0.5)	2211
A10 (opacious A1)	2007	A10 (opacious A1)	2212
A20 (opacious A2)	2008	A20 (opacious A2)	2213
Inc (Incisal)	2010	Inc (Incisal)	2215
BW (Bleaching White)	2011	BW (Bleaching White)	2216

F03 – moderate flow (Low Flow)

Shade	PN 2.2 g	Shade	PN 0.21 g
SYRINGE		TIPS, 20 pcs.	
A1	2014	A1	2219
A2	2015	A2	2220
A3	2016	A3	2221
A3.5	2017	A3.5	2222
A4	2018	A4	2223
B1	2066	B1	2224
B2	2067	B2	2225
C2	2068	C2	2226
D2	2069	D2	2227
SA0.5	2013	SA0.5	2218
A20 (opacious A2)	2020	A20 (opacious A2)	2229
A30 (opacious A3)	2021	A30 (opacious A3)	2230
Inc (Incisal)	2023	Inc (Incisal)	2232
BW (Bleaching White)	2024	BW (Bleaching White)	2233
MI (Milky)	2026	MI (Milky)	2235
CV (Cervical)	2027	CV (Cervical)	2236

Standard Kit



Standard Kit

- 1 syringe each at 2.2 g Beautifil Flow Plus F00 and F03, shade A2
- 1 syringe each at 2.2 g Beautifil Flow Plus F00 and F03, shade A3

PN	Quantity	Product
2000S	1 pc.	Beautifil Flow Plus Standard Kit
2000-2027, 2062-2069	2.2 g	Beautifil Flow Plus - syringe
2200-2236	0.21 g	Beautifil Flow Plus - tips, 20 pcs.