

# Bur Stands

For systematic storage and preparation in the dental practice

The Bur Block and the Bur Station II are ideal bur stands for the storage of contra angle (CA) and friction grip (FG) instruments. The organisers are handy and aid easy transportation of up to 12 FG and 6 CA instruments. Depending on the indication, various individual sets can be assembled.



## Bur Block

The lockable lid protects the instruments against spillage and contamination, the open sides ensure total sterilisation and disinfection of up to 18 instruments (12 FG and 6 CA). Depending on the procedure (e.g. cavity preparation), individual colour coded instrument sets can be assembled.

- Aluminium housing with lockable lid
- 5 colours: black, blue, gold, green, red
- Dimension: 4.4 x 3.1 x 2.6 cm
- For 12 FG and 6 CA instruments



## Bur Station II

The Bur Station II for storage and sterilisation of up to 17 instruments (11x FG and 6x CA) is made from multi functional resin, and offers a high thermal resistance (< 140 °C), thus being very durable and resistant to chemicals. It is available in two colours.

- Multi-functional resin stand with lockable lid
- 2 colours: yellow, light blue
- Dimension: 4.2 x 2.2 x 3 cm
- For 11 FG and 6 CA instruments

## BENEFITS

Can be closed securely

Can be disinfected and sterilised in an autoclave

PN	Quantity	Product
92BK	1 pc.	Bur Block, black
92BL	1 pc.	Bur Block, blue
92GL	1 pc.	Bur Block, gold
92GN	1 pc.	Bur Block, green
92MN	1 pc.	Bur Block, red
0796	1 pc.	Bur Station II, yellow
0797	1 pc.	Bur Station II, light blue