

CompoMaster Finishing & Polishing Kit

For systematic finishing and polishing of composite restorations

The kit contains all the necessary rotary instruments for systematic finishing and polishing of composite restorations in the dental practice.

Contents

T&F Hybrid Points: 1 each 7106, 7406, 7901, LT1
 CompoMaster Coarse: 1 each knife, cup, minipoint
 CompoMaster: 1 each knife, cup, minipoint

PN	Quantity	Product
0331	1 pc.	CompoMaster Finishing & Polishing Kit, CA-FG shank combination

