

Composite Polishing Kit



For finishing and polishing of composite restorations

The kit contains all instruments for contouring, finishing and polishing all types of composite restorations: Dura-Green stones to reduce excess material, Dura-White stones for finishing, CompoSite polishers for prepolishing and CompoSite Fine instruments for high gloss without damage to the enamel.

Contents

Dura-Green: 1 each CN1, FL2

Dura-White: 1 each CN1, FL2

CompoSite: 2 each bullet, knife, wheel

CompoSite Fine: 1 each cup, disk

PN	Quantity	Product
0310	1 pc.	Composite Polishing Kit, CA shank

