

# Dura-Polish ■ Dura-Polish DIA

## Polishing paste system

Dura-Polish is impregnated with up to 73 % by weight with aluminum oxide particles, and efficiently achieves a very homogeneous surface when prepolishing composite, alloys, denture acrylics and PMMA, especially in difficult to reach areas.

Dura-Polish DIA is filled to an extremely high level (67 % by weight) with industrial diamond particles. Due to the ultra fine particle size of  $< 1 \mu\text{m}$  this polishing paste achieves an excellent and consistent high gloss with composite and porcelain.



### Indications

#### ■ DURA-POLISH

Prepolishing and polishing of composite, alloy and denture acrylics

#### ■ DURA-POLISH DIA

High-lustre polishing of composite and porcelain, zirconium and alloys



Min<sup>-1</sup> Max. 10,000

Not suitable for intraoral use!

## BENEFITS

Adhesive character of paste ensures maximum surface contact

Viscosity of paste minimises splattering

Economical to use

Excellent polishability

PN	Quantity	Product
0552	1 jar	Dura-Polish, 20 g
0554	1 jar	Dura-Polish DIA, 5 g

