

# Gingi Curettage Diamonds

For crown and bridge preparation and grinding of tooth stumps

The gold coloured standard grit FG diamonds are designed to establish a chamfered finish to the tooth surface in crown and bridge preparations, giving an ideal profile for precise impression taking. They are made with a 3° axial side taper, a 45° tip angle to the axis and a 1.52 mm point radius.

PN	C061-1	C062-1	C063-1	C081-1	C082-1	C083-1	C084-1	C085-1	C101-1	C102-1	C103-1
ISO No.	012	014	017	013	016	018	019	021	014	015	021
Head length (mm)	6.0	6.0	6.0	8.0	8.0	8.0	8.0	8.0	10.0	10.0	10.0



Min<sup>-1</sup> Recom. 160,000 - 450,000 / CP 0 - 0.8 N

PN	Quantity	Product
C061-1 - C103-1	6 pcs.	Gingi Curettage Diamonds, FG shank
C200	1 pc.	Gingi Curettage Diamonds Kit

## Indications

CHAIRSIDE

Crown and bridge preparation, grinding of tooth stumps

## Range of Application

Suitable for the following materials:  
Tooth enamel and dentine



## BENEFITS

- Exact, even distribution of the diamond particles
- Precise concentricity
- Low vibration performance
- High cutting performance with low pressure and secure cooling

## info

The colour coding of the order number indicates the individual maximum speed:

- Cyan max. 450,000 rpm
- Blue max. 300,000 rpm
- Orange max. 160,000 rpm



The Gingi Curettage Diamonds Kit contains the six most popular shapes C062-1, C081-1, C083-1, C084-1, C102-1, C103-1.