



# Glaslonomer FX-II

All-purpose glass ionomer, specially created for the minimally invasive dentistry

Glaslonomer FX-II is an universal and modellable glass ionomer for general dentistry. It was specially created for minimally invasive treatments.



## Indications

- Minimally invasive treatments
- Final restorative for deciduous teeth
- Geriatric restorative for class I, II, III, V cavities and cervical erosions
- Final restorative for class I, II of adult dentition in non-stress bearing areas
- Temporary restorative for heavy stressed cavities
- Core build-up and liner
- Pit and fissure sealing

PN	Quantity	Product
3192-3195	1 pc.	<b>Glaslonomer FX-II Cement Set</b> powder 15 g, liquid 8 ml, cocoa butter 1 g, matrix strips, mixing pad, spatula, A2: 3192, A3: 3193, A3.5: 3194, B2: 3195
3203	1 pc.	<b>Glaslonomer FX-II Cement Mini-Set</b> shade A3, powder 6 g, liquid 2.8 ml, spoon