



HY-Bond Polycarboxylate Cement

Luting cement

HY-Bond Polycarboxylate Cement contains a unique tannin-fluoride additive to help sealing the dentinal tubules.



Indications

- Direct luting of crowns, bridges and inlays
- Lining of cavities and lesions

PN	Quantity	Product
1160	1 pc.	HY-Bond Polycarboxylate Cement Set powder 60 g, liquid 33 ml, spoon
1161	60 g	HY-Bond Polycarboxylate powder
1162	33 ml	HY-Bond Polycarboxylate liquid

BENEFITS

- Chemical bond of tooth structure and metal
- High adhesive strength
- Radiopaque
- Contains tannin-fluoride