

# HY-Bond Temporary Cement – Soft



## Temporary polycarboxylate cement

With this temporary cement type the polycarboxylate cement base was combined with a tannin-fluoride additive. It is ideally suited for provisional bonding.



### Indications

- Temporary cementing of provisionals
- Indirect pulp capping, temporary sealing after cavity preparation
- Temporary cementing of deciduous cavities

## BENEFITS

- Short setting time
- Eugenol-free
- Odourless and tasteless
- Easy to remove
- Easy to mix and convenient to handle
- Contains a patented tannin-fluoride additive (HY-Agent)

PN	Quantity	Product
1173	1 pc.	HY-Bond Temporary Cement - Soft Set powder 125 g, liquid 60 ml, spoon