

LAMINA VEST II

Phosphate bonded investment for minimally invasive porcelain techniques



Lamina Vest II is a phosphate bonded investment for the production of refractory dies. Developed for the fabrication of porcelain inlays, veneers and crowns in the sintering technique, Lamina Vest II offers high strength and edge stability.

BENEFITS

Suitable for Vintage PRO, Vintage MP, Vintage Halo and Vintage ZR

Excellent fit of the porcelain restorations

Smooth surfaces

High strength and edge stability

High dimensional stability, even after repeated firings

Easy divesting

Indications

For all porcelains with a CTE range of 9.5 to $12.5 \times 10^{-6}K^{-1}$ ($25\text{ }^{\circ}C$ - $500\text{ }^{\circ}C$)

| PN | Quantity | Product |
|------|----------|--|
| 6977 | 1.2 kg | Lamina Vest II powder (40 x 30 g bags) |
| 6978 | 200 ml | Lamina Vest II liquid (with measuring cup) |
| 6972 | 2 pcs. | Refractory marker |

