



# Mandrels

## Mandrels for unmounted wheels and disks

The mandrels from stainless steel are available with hand piece, contra-angle and friction grip shank (CA, FG). Single use resin mandrels are available in contra-angle. Further information regarding handling can be found on page 6.

PN	Quantity	Product
0435	6 pcs.	CA mandrel for OneGloss polishers
0439	6 pcs.	CA mandrel for Super-Snap disks, Super-Snap SuperBuff, Super-Snap Buff Disks
0438	6 pcs.	FG mandrel for Super-Snap disks, Super-Snap SuperBuff, Super-Snap Buff Disks
0440	100 pcs.	CA mandrel, resin, for Super-Snap disks



0435

0439

0438

0440