

Mandrels

Mandrels for unmounted wheels and disks

The mandrels from stainless steel are available with hand piece, contra angle and friction grip shank (HP, CA, FG). Single use resin mandrels are available in contra angle. Further information regarding handling can be found on page 6.

PN	Quantity	Product
5011	12 pcs.	HP mandrel for all unmounted wheels and disks
0435	6 pcs.	CA mandrel for OneGloss polishers
0439	6 pcs.	CA mandrel for Super-Snap disks, Super-Snap SuperBuff, Super-Snap Buff Disks
0438	6 pcs.	FG mandrel for Super-Snap disks, Super-Snap SuperBuff, Super-Snap Buff Disks
0440	100 pcs.	CA mandrel, resin, for Super-Snap disks



5011



0435



0439



0438



0440