

Veracia

Composite teeth for full and partial dentures

The sculptured tooth shape and anatomical design of Veracia anterior and posterior teeth guarantee multifunctional purposes and provide an unlimited suitability for all standard setup techniques of full and partial dentures.

Veracia artificial teeth cannot be compared to conventional acrylic teeth. Based on the technology of the proven Solidex c & b composite system, the teeth consist of micro filled hybrid composite. By combining organic and inorganic fillers in the composite matrix, an extraordinary abrasion resistance is achieved.



Indications

- Full denture prosthetics
- Implant restorations
- Telescopic restorations
- Attachment work
- Cr-Co dentures

BENEFITS

Lively surface structure, matching the morphology of natural teeth

Multifunctional use for all standard set-up techniques

Limited variety of moulds for easier choice of anterior and posterior teeth as well as economical storage

Excellent match with the VITA* Classical Shade System

Setup techniques

- Tooth to tooth occlusion
- Tooth to two teeth occlusion
- Balanced and lingualised occlusion
- Incisal and canine-guided occlusion
- Sequential guidance

Veracia - Accessories

PN	Quantity	Product
4061	1 pc.	Dental Measure II

Veracia - Cabinet

Contents

576 anteriors and 384 posteriors of all common shapes in the VITA* Classical shades A2, A3, A3.5, B2 and B3

PN	Quantity	Product
AOC001	1 pc.	Veracia Cabinet





Veracia - Living Mould Chart

Contents

15 upper and 8 lower anterior tooth shapes (138 anteriors) and 3 upper and lower posteriors (48 posterior teeth) in the shade A3

PN	Quantity	Product
AOLM01	1 pc.	Veracia Living Mould Chart



Veracia - 28 Teeth Set

Contents

3 different types of anterior shapes in combination with the posterior shape "M" in the VITA* Classical shade A3

PN	Anterior upper / lower jaw	Posterior upper / lower jaw
A3STA	A3ST5 / A3A5	A3M0 / A3M1
A3LSL	A3LS5 / A3LA5	A3M0 / A3M1
A3LTL	A3LT5 / A3LA5	A3M0 / A3M1

Veracia - Tooth Cards

All sets are available in the VITA* Classical shades A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4 and in the SHOFU whitening shade W1. Please ask for our mould chart and order information.

6 pcs. each per tooth card



Veracia - Upper Anterior Teeth

Available in the shapes O4, O5, O6, LO5, LO6, S4, S5, S6, LS5, LS6, ST4, ST5, ST6, LT5, LT6



Veracia - Lower Anterior Teeth

Available in the shapes NA4, NA5, NA6, A4, A5, A6, LA5, LA6



8 pcs. each per tooth card

Veracia - Posterior Teeth

Available in the shapes S upper, S lower, M upper, M lower, L upper, L lower

