

Veracia



Fully anatomical composite teeth for full and partial dentures

Veracia anterior and posterior composite teeth are manufactured using state-of-the-art CAD/CAM technology, based on the anatomical parameters of natural teeth – for exceptional, impressive aesthetics and lifelike surface structures.

The teeth are multifunctional, can be used in all standard setup techniques for complete and partial dentures without limitations, and meet the ever-increasing demand of patients for aesthetic and functional prosthodontic solutions.



Indications

- Full denture prosthetics
- Implant restorations
- Telescopic restorations
- Attachment work
- Cr-Co dentures

BENEFITS

- Lively surface structure, matching the morphology of natural teeth
- Multifunctional use for all standard set-up techniques
- Limited variety of moulds for easier choice of anterior and posterior teeth as well as economical storage
- Excellent match with the VITA® Classical Shade System

Setup techniques

- Tooth to tooth occlusion
- Tooth to two teeth occlusion
- Balanced and lingualised occlusion
- Incisal and canine-guided occlusion
- Sequential guidance

Veracia - Accessories

PN	Quantity	Product
4061	1 pc.	Dental Measure II

Veracia - Cabinet

Contents

576 anteriors and 384 posteriors of all common shapes in the VITA® Classical shades A2, A3, A3.5, B2 and B3

PN	Quantity	Product
AOC001	1 pc.	Veracia Cabinet





Veracia - Living Mould Chart

Contents

15 upper and 8 lower anterior tooth shapes (138 anteriors) and 3 upper and lower posteriors (48 posterior teeth) in the shade A3

PN	Quantity	Product
AOLM01	1 pc.	Veracia Living Mould Chart



Veracia - 28 Teeth Set

Contents

3 different types of anterior shapes in combination with the posterior shape "M" in the VITA* Classical shade A3

PN	Anterior upper / lower jaw	Posterior upper / lower jaw
A3STA	A3ST5 / A3A5	A3M0 / A3M1
A3LSL	A3LS5 / A3LA5	A3M0 / A3M1
A3LTL	A3LT5 / A3LA5	A3M0 / A3M1

Veracia - Tooth Cards

All sets are available in the VITA* Classical shades A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4 and in the SHOFU whitening shade W1. Please ask for our mould chart and order information.

6 pcs. each per tooth card




Veracia - Upper Anterior Teeth

Available in the shapes O4, O5, O6, LO5, LO6, S4, S5, S6, LS5, LS6, ST4, ST5, ST6, LT5, LT6



Veracia - Lower Anterior Teeth

Available in the shapes NA4, NA5, NA6, A4, A5, A6, LA5, LA6



8 pcs. each per tooth card

Veracia - Posterior Teeth

Available in the shapes S upper, S lower, M upper, M lower, L upper, L lower

