



# VINTAGE LD Press

Ingots for the production of all-ceramic restorations

Vintage LD Press is a lithium disilicate glass ceramic in form of ingots for pressing inlays, onlays, veneers, crowns and anterior bridges. The press ingots have a high mechanical strength and are available in four different translucencies.

## ■ T (Trans)

high translucency in three shades  
for producing smaller restorations like inlays,  
onlays and veneers in the staining technique

## ■ MT (Middle Trans)

middle translucency in all VITA\* Classical shades  
for the cut-back, layering or staining technique

## ■ LO (Low Opacity)

low opacity in five colours  
for individual abutments or slightly discoloured  
preparations

## ■ MO (Middle Opacity)

middle opacity in three shades  
for production of frameworks for strongly  
discoloured preparations



## BENEFITS

Long-term stability due to high  
mechanical strength

Simple and efficient handling

Four different translucencies for  
natural aesthetics, independent of the  
die colour

Can be combined with Vintage LD  
layering porcelain and Vintage Art LF  
stains



## Indications

For pressing crowns, inlays, onlays, veneers and  
anterior bridges (up to three units without molars)

## VINTAGE LD Press - Sets and Single Components

Basic Set



### Basic Set

For pressing metal-free lithium disilicate restorations in the shades A2, A3 and B2 as well as ingots of different opacities for the layering and staining technique

#### Vintage LD Press

- 40 ingots, 3 g each: 5 pcs. each in the shades MT-A1, MT-A2, MT-A3, MT-B2, L0-1, L0-2, T-1, T-2

#### Ceravety Press & Cast

- 5x powder, 100 g each • liquid, 100 ml

### Trial Kit

5 ingots of different shades and translucencies for the introduction of the Vintage LD press technique

- 5 ingots, 3 g each: 2x MT-A2, 2x MT-A3, 1x T-1



PN	Quantity	Product
8501	1 pc.	Vintage LD Press Basic Set
8542	1 pc.	Vintage LD Press Trial Kit
8510-8540	1 pc.	Vintage LD Press ingots, 5 pcs.

Shade	PN 5 pcs.
VINTAGE LD PRESS	
MT-W1	8527
MT-W2	8528
MT-W3	8529
MT-A1	8510
MT-A2	8511
MT-A3	8512
MT-A3.5	8513
MT-A4	8514
MT-B1	8515
MT-B2	8516
MT-B3	8517
MT-B4	8518
MT-C1	8519
MT-C2	8520
MT-C3	8522

Shade	PN. 5 pcs.
VINTAGE LD PRESS	
MT-C4	8523
MT-D2	8524
MT-D3	8525
MT-D4	8526
T-1	8530
T-2	8531
T-3	8532
L01	8533
L02	8534
L03	8535
L04	8536
L05	8537
M01	8538
M02	8539
M03	8540

