

# Dressing Tools



For shaping and cleaning of abrasives and polishers

For the precise work with abrasives and polishers shapes often need to be customised. Besides, sharp edges or rounded parts may be worn, so that the instrument is no longer used, even though there still is enough material on the shank.

The diamond impregnated Abrasives Dresser for ceramic-bonded stones and polishers and the SHOFU Adjuster made of alumina for the Dura-Green DIA instruments enable easy shaping of edges and rounded parts, as well as cleaning of the instruments.



PN	Quantity	Product
0781	1 pc.	Abrasives Dresser, 100 x 15 x 4 mm
0504	1 pc.	SHOFU Adjuster, 85 x 17 x 10 mm