

# AcryPoint-E

For PMMA resins



The AcryPoint-E silicone polishers are designed for adjusting, finishing and polishing of acrylic dentures, all acrylic appliances and acrylic impression trays without the need to use additional polishing pastes.



## Indications

### LABORATORY

- The coarse dark grey AcryPoint-E are used for adjustments.
- The brown AcryPoint-E with medium grit are designed for finishing and prepolishing.
- The light grey AcryPoint-E with fine grit are recommended for polishing.



## highlights

- Three well-matched polishing steps
- Optimal surface polishing of PMMA
- Optimal contouring and prepolishing of acrylic bite splints in one step
- Exact match of abrasive grains and binder
- High gloss polishing without the need for an additional polishing agent
- Low vibration performance
- Low heat development



	Shape	BP1	PC2
	ISO No.	100	050
Dark grey	HP 6 pcs.	SE0009	
	HP 12 pcs.		SE0006
Brown	HP 6 pcs.	SE0010	
	HP 12 pcs.		SE0007
Light grey	HP 6 pcs.	SE0011	
	HP 12 pcs.		SE0008



**Min<sup>-1</sup>** Recom. 5,000- 10,000 / Max. 15,000 / CP 3-5 N  
Do not use intraorally!

PN	Quantity	Product
SE0006 - SE0008	12 pcs.	AcryPoint-E, shape PC2, HP shank
SE0009 - SE0011	6 pcs.	AcryPoint-E, shape BP1, HP shank