

All Ceramic Preparation Kit



Preparation system for all-ceramic restorations

This assortment of 17 instruments is ideally suited for the preparation of anterior teeth, premolars and molars, as well as inlay and veneer preparations. The shapes of the abrasives are matched to the different preparation designs, ensuring optimal congruence of shapes.

Contents

Robot Points - 1 each:

808V-1, 0825-1, 825V-1, 0826-1, 826V-1, 0833-1,
833V-1, 0836-1, 836V-1, 0839-1, 839V-1, 0893-1,
893V-1, 0896-1, 0897-1, 0899-1, 899V-1

PN	Quantity	Product
0955	1 pc.	All Ceramic Preparation Kit

