

Amalgam Polishing Kit



Polishing kits designed for restorations made of amalgam or precious alloys, also suitable for all composites and ceramics

The kits are available in contra-angle or friction grip shanks and contain the instruments needed for all polishing steps up to high gloss.

Contents

Cup: 3x Brownie, 2x Greenie, 1x Supergreenie

Minipoint: 3x Brownie, 2x Greenie, 1x Supergreenie

PN	Quantity	Product
0308	1 pc.	Amalgam Polishing Kit, CA shank
0309	1 pc.	Amalgam Polishing Kit, FG shank

