

CeraMaster Coarse ■ CeraMaster



For ceramic restorations

These silicone polishers are integrated with pure, ultrafine diamond particles. They are characterised by a long service life and designed for quick finishing and polishing of ceramic restorations, to provide a smooth surface. The combination of the different grits of CeraMaster Coarse (dark-grey polishers) and the finer CeraMaster polishers (blue and white banded shank) is an efficient concept for homogeneous surface treatment.




	Shape	Bullet	Knife	Cup	Minipoint	KN7	WH6	Unmounted wheel
	ISO No.	040	090	065	030	125	120	220
CeraMaster Coarse	HP	126C			129C	127C	128C	132C
	CA	121C	122C	123C	124C			
CeraMaster	HP	0126			0129	0127	0128	
	CA	0121	0122	0123	0124			


Indications

■ CeraMaster Coarse polishers are designed for a quick finishing and prepolishing of ceramic restorations in one step only. They work excellently even when adjusting margins of all-ceramic frame-works made of alumina or zirconia.

■ For super-polishing CeraMaster polishers are recommended. When using the CeraMaster polishers there is no need to use an additional polishing paste, as the high shine is produced by using the diamond-impregnated silicone polishers themselves.

PN	Quantity	Product
INDIVIDUAL ARTICLES		
121C-124C, 126C-129C	3 pcs.	CeraMaster Coarse, HP and CA shank
132C	1 pc.	CeraMaster Coarse, unmounted wheel + HP mandrel
0121-0124, 0126-0129	3 pcs.	CeraMaster, HP and CA shank
ASSORTMENTS		
130C	1 box	CeraMaster Coarse, HP shank, 1 each bullet, minipoint, KN7, WH6
125C	1 box	CeraMaster Coarse, CA shank, 1 each bullet, knife, cup, minipoint
0130	1 box	CeraMaster, HP shank, 1 each bullet, minipoint, KN7, WH6
0125	1 box	CeraMaster, CA shank, 1 each bullet, knife, cup, minipoint

 **Polishers:**
Min⁻¹ Recom. 10,000 - 15,000 / Max. 20,000 / CP 1.5 N

 **Unmounted wheel:**
Min⁻¹ Recom. 5,000 - 10,000 / Max. 15,000 / CP 1.5 N

