CompoSite - CompoSite Fine





For composite

CompoSite instruments are designed for a perfect finish of composite restorations without scratching the surrounding enamel. CompoSite polishers are used for finishing and prepolishing, CompoSite Fine (with white banded shank) for polishing. The flexible and adaptable cups and disks will polish subgingival and interproximal surfaces safely while pear points and midipoints are recommended for the polish of occlusal and lingual surfaces.





highlights

Silicone polishers with integrated abrasive grains

Matched hardness of silicone

Optimal ratio of speed and contact pressure

No damage to enamel

No polishing paste is needed

Long-lasting tool life



CA shank:

Min-1 Recom. 10,000-12,000 / Max. 20,000 / CP 0.5 N



FG shank:

Min⁻¹ Recom. 10,000-12,000 / Max. 15,000 / CP 0.5 N

Indications

DENTAL PRACTICE

For finishing and polishing of micro and hybrid composite restorations such as inlays, onlays, crowns and bridges in the patient's mouth as well as for removing cements used for orthodontic brackets



0191-0299, 196B, 296B 12 pcs. CompoSite / CompoSite Fine, CA and FG shank

190b, 290b CA and FG Sharik



	Shape	Bullet	Knife	Wheel	Disk	Cup	Pear-Point	IC	Midipoint
	ISO No.	040	090	090	160	080	033	090	040
CompoSite	CA	0191	0192	0193	0194	0195	0196	0197	0198
	FG						196B		0199
CompoSite Fine	CA	0291	0292	0293	0294	0295	0296	0297	0298
	FG						296B		0299