

Composite Technique Kit



For composite restorations

Quickly, easily and without the need for additional materials, CompoSite and CompoSite Fine polish the entire surface structure of a white filling in only two steps. Simultaneously the surface is sealed to prevent bacterial infiltration.

Contents

- Dura-White: 1 each CN1, FL2
- CompoSite: 1 each bullet, knife, wheel, disk, cup, midpoint, pear-point, IC
- CompoSite Fine: 1 each bullet, knife, wheel, disk, cup, midpoint, pear-point, IC

PN	Quantity	Product
0322	1 pc.	Composite Technique Kit, CA shank

