

# Dura-Polish ■ Dura-Polish DIA



Polishing paste system for composite and porcelain

Dura-Polish is impregnated with fine alumina particles (more than 73 % by weight) and efficiently achieves a very homogeneous surface when prepolishing composite, alloys, denture acrylics and PMMA, especially in difficult to reach areas.

Dura-Polish DIA is filled to an extremely high level (67 % by weight) with industrial diamond particles. Due to the ultrafine particle size of < 1 µm this paste polishes composites and porcelain to an excellent and durable gloss in a minimum of time.



## highlights

Good adhesion, reliably sticks to the surface

Firm consistency, for minimum risk of splattering

Economical to use

Excellent polishing properties

PN	Quantity	Product
0552	1 jar	Dura-Polish, 20 g
0554	1 jar	Dura-Polish DIA, 5 g

### Indications

#### ■ DURA-POLISH

Prepolishing and polishing of composite, alloy and denture acrylics

#### ■ DURA-POLISH DIA

High-lustre polishing of composite and porcelain, zirconia and alloys



Min<sup>-1</sup> Max. 10,000  
Not suitable for intraoral use!