

Enamel Adjustment Kit



For adjustment and systematic polishing of tooth enamel and porcelain

The Enamel Adjustment Kit is indispensable to dental practice for very fine adjustments of natural dentition or porcelain restorations and subsequent super-polish in the patient's mouth. The choice of shapes is particularly well suited for working on difficult to reach areas, such as interproximal areas or occlusal surfaces.

Contents

Dura-White: 1 each FL2, RD1, TC1
Ceramisté Standard: 1 each knife, wheel, bullet
Ceramisté Ultra: 1 each knife, wheel, bullet
Ceramisté Ultra II: 1 each knife, wheel, bullet



PN	Quantity	Product
0307	1 pc.	Enamel Adjustment Kit, CA shank