

# Merssage ■ Pressage



## Prophylaxis pastes for all indications

These prophylaxis pastes are available in different types. The fluoride containing Merssage pastes remove plaque and provide a smooth and polished tooth surface. The fluoride-free Pressage paste is designed for cleaning tooth surfaces before bleaching or pit & fissure sealing.

### Indications

#### ■ MERSSAGE REGULAR

RDA 170-180, with fluoride

Merssage Regular – with peppermint taste – is used for the removal of discolouration and at the beginning of the prophylaxis treatment, as well as for polishing teeth surfaces after the removal of orthodontic appliances.

#### ■ MERSSAGE FINE

RDA 40-50, with fluoride

Merssage Fine – with lemon taste – is used for the removal of plaque and extrinsic stains on the tooth surfaces and for polishing after scaling or bleaching.

#### ■ MERSSAGE PLUS

RDA < 10, with fluoride


Merssage Plus – with orange taste – is used for polishing after scaling or bleaching, for prevention of reattached stains as well as for cleaning of restorative materials and polishing the surface.

#### ■ PRESSAGE

RDA 170-180, without fluoride

The tasteless Pressage paste is used for the surface cleaning before bleaching or application of etchants, bonding agents or fissure sealants. It does not contain fluoride, as this can weaken the bonding strength and the bleaching effect.



 Min<sup>-1</sup> Max. 2,000  
Do not operate with water spray!

PN	Quantity	Product
0572	1 tube	Merssage Regular, 40 g
0573	1 tube	Merssage Fine, 40 g
0575	1 tube	Merssage Plus, 38 g
0574	1 tube	Pressage, 40 g