

OneGloss M ■ OneGloss ■ OneGloss PS



One-step finishers and polishers for composites and enamel

These one-step silicone polishers, impregnated with alumina, allow you to perfectly finish and polish all types of composite restorations, simply by varying the contact pressure step by step from finishing to polishing – without instrument change and polishing paste!

The product line includes four shapes which are optimally matched to the surfaces to be treated. The OneGloss M instruments, with heads mounted on stainless-steel contra-angle shanks, are sterilisable. The unmounted OneGloss instruments and the OneGloss PS instruments, with heads mounted on sturdy plastic mandrels, are disposable.

Indications

■ DENTAL PRACTICE

Finishing and polishing of all composite restorations and tooth enamel in one step, removal of excess cement without marring delicate inlays, removal of surface stains and polishing of tooth enamel after professional cleaning



highlights

- Silicone polishers with alumina
- Optimal ratio of speed and contact pressure
- No need to change instruments
- No polishing paste is necessary
- Optimal abrasive performance
- Low vibration action

OneGloss M / OneGloss:
Min⁻¹ Recom. 3,000 - 10,000 / Max. 15,000 /
 CP: 0.3 N for polishing / 1.0 N for finishing

OneGloss PS:
Min⁻¹ Recom. 3,000 - 5,000 / Max. 5,000 /
 CP: 0.3 N for polishing / 1.0 N for finishing



Shape	Cup	Midpoint	IC	Minipoint	Mandrel
ISO No.	060	060	100	030	
OneGloss M	0187	0188	0189	0186	
OneGloss	0181	0182	0183		0435
OneGloss PS	0176	0177	0178		

info

Use intermittent water spray when applying high pressure or polishing enamel.

OneGloss M polishers are sterilisable.

OneGloss and OneGloss PS polishers are disposable.

PN	Quantity	Product
OneGloss M		
0186-0189	12 pcs.	OneGloss M, mounted, CA shank shapes: cup, midpoint, IC, minipoint
OneGloss		
0180	1 pc.	OneGloss Kit 60 unmounted polishers (20 each cup, midpoint, IC) and 3 CA mandrels
0181-0183	1 pc.	OneGloss Refill 50 pcs., unmounted, shapes: cup, midpoint, IC
0435	6 pcs.	CA mandrel
OneGloss PS		
0175	1 pc.	OneGloss PS Kit 50 pcs.: 15x cup, 15x IC, 20x midpoint, all individually wrapped
0176-0178	50 pcs.	OneGloss PS, all individually wrapped, shapes: cup, midpoint, IC