



Post Stand P ■ Post Set

Tray and pin set

Post Stand P

Two crossing central ditches extending down to the bottom of the Post Stand and more than 100 holes allow – combined with the matching pins – an exceptionally fast and simple positioning of the restorations, from single crowns up to 14 unit bridges. Dimensions: ø 60 mm / h 23.5 mm



Post Set

The set includes five different shapes, 3 each, for any kind of restoration:

- **Anterior Post A** (cone shaped) and **Anterior Post B** (wedge shaped) for anterior crowns
- **Posterior Post A** for premolars
- **Posterior Post B** for molars
- **Free Post** can be grinded to any desired shape

By decentralised shanks and rotatability the space between the posts can be individually adjusted.



Indications

For stable placement of restorations in dental lab light-curing units

PN	Quantity	Product
5059	1 pc.	Post Stand P (tray)
5058	1 pc.	Post Set (pins)